



KONICA MINOLTA

Inspection & Calibration system for CRI Illuminance meter CL-70F

KONICA MINOLTA, INC.
Instrument Service Dept.
Sensing Business Unit

		SECONIC CORPORATION	KONICA MINOLTA Primary authorized service facility (Inspection only)
National Institution		National Institute of Advanced Industrial Science and Technology (Japan)	
International standard		↓	
Designated calibration laboratories or JEMIC	National primary standard:	Japan Electric Meters Inspection Corporation	
Accredited laboratories	National secondary standard:	↓	
		JCSS-0039 Japan Electric Meters Inspection Corporation	
Factory	In-house standard	↓	
		Distribution Temperature / Luminous Intensity Standard Lamp	
		↓	
		Illuminance Distribution Temperature	Master body
		↓	↓
		Illuminance Distribution Temperature	Light Source
		↓	↓
Product		CRI Illuminance meter CL-70F	CRI Illuminance meter CL-70F

* This is applicable at the calibration by service facilities of KONICA MINOLTA authorization

JCSS-0039: Japan Calibration Service System is the calibration laboratory accreditation system in Japan. 0039 is the registration number of JEMIC.